


A Rotating Grater-like Cutting Tool to Remove Hard Cemented Deposits in Heart Blood Vessels without Damaging Soft Vessel Walls

Tsuchiya Laboratory, Institute of Industrial Science, The University of Tokyo



Background




PTCRA (percutaneous transluminal coronary rotational atherectomy)



Design of rotating grater-like cutting tool

- FR1: Remove hard cemented deposits, i.e. calcified atheroma.
 - FR2: Do not damage the soft vessel walls.
 - FR3: Do not let tool tips wear or fall off.
 - FR4: Finish cutting the calcified atheroma within 30 seconds or less.
 - FR5: Be visible with X-ray.
 - FR6: Stop the cutting upon excessive cutting force.
- (FR: Functional Requirement)

Fabrication of the grater tool



Bullet shaped tool (Steel)



Type	#	material	height [m]	pitch [mm]	orientation
Grater height	1	anodized aluminum	30	0.2	dent first
	2		20	0.2	dent first
	3		20	0.2	bulge first
	4		5	0.2	bulge first
	5		5	0.1	bulge first
	6	steel*1	20	0.2	bulge first
	7	brass*2	20	0.2	dent first
	8	brass	-20	0.2	dent only
	9	brass	0	-	flat
Diamond height	10	nickel plated brass	30	-	-
	11		20	-	-
	12		15	-	-
	13		10	-	used in surg.
	14		7	-	-

*1 quenched tool steel *2 nickel plated brass


Table 1: Used tools in the experiments.

Experiments

Experimental set up



The tool in rotating is pressed down to 500µm against the egg underwater



(a) w/ grater tool (b) w/ diamond tool

Best tools for cutting eggshell

Type	#	material	height [m]	pitch [mm]	orientation	#	removal rate of eggshell cutting (FR1&FR5)	scratch on cut egg white (FR2)	used tool surface aft. 3min. eggshell cutting (FR3)
Grater height	1	anodized aluminum	30	0.2	dent first	1	0.095 [mm³/s]	visible	no wear
	2		20	0.2	dent first	2	0.089	visible	no wear
	3		20	0.2	bulge first	3	0.074	not visible	no wear
	4		5	0.2	bulge first	4	0.010	not visible	no wear
	5		5	0.1	bulge first	5	0.025	slightly visible	no wear
	6	steel*1	20	0.2	bulge first	6	0.069	not visible	local wear on bulges
	7	brass*2	20	0.2	dent first	7	0.008	visible	hard wear on bulges
	8	brass	-20	0.2	dent only	8	0	slightly visible	-
	9	brass	0	-	flat	9	0	not visible	-
Diamond height	10	nickel plated brass	30	-	-	10	0.15	visible	-
	11		20	-	-	11	0.13	visible	-
	12		15	-	-	12	0.10	visible	-
	13		10	-	used in surg.	13	0.087	not visible	particles dropped
	14		7	-	-	14	0.072	not visible	particles dropped

*1 quenched tool steel *2 nickel plated brass

Table 1: Used tools in the experiments.

Table 2: Experimental result (screening columns meet the FRs).

Grater : 20µm high is the best
Diamond: 7µm high is the best

Water flow influence

Type	#	material	height [m]	pitch [mm]	orientation	#	removal rate of eggshell cutting (FR1&FR5)	scratch on cut egg white (FR2)	used tool surface aft. 3min. eggshell cutting (FR3)
Grater height	1	anodized aluminum	30	0.2	dent first	1	0.095 [mm³/s]	visible	no wear
	2		20	0.2	dent first	2	0.089	visible	no wear
	3		20	0.2	bulge first	3	0.074	not visible	no wear
	4		5	0.2	bulge first	4	0.010	not visible	no wear
	5		5	0.1	bulge first	5	0.025	slightly visible	no wear
	6	steel*1	20	0.2	bulge first	6	0.069	not visible	local wear on bulges
	7	brass*2	20	0.2	dent first	7	0.008	visible	hard wear on bulges
	8	brass	-20	0.2	dent only	8	0	slightly visible	-
	9	brass	0	-	flat	9	0	not visible	-
Diamond height	10	nickel plated brass	30	-	-	10	0.15	visible	-
	11		20	-	-	11	0.13	visible	-
	12		15	-	-	12	0.10	visible	-
	13		10	-	used in surg.	13	0.087	not visible	particles dropped
	14		7	-	-	14	0.072	not visible	particles dropped

Water flow doesn't cut the shell.

not visible at 200,000 rpm → visible at 7,000 rpm


Water flow prevents the egg white scratch.

Wear of the tool

Type	#	material	height [m]	pitch [mm]	orientation	#	removal rate of eggshell cutting (FR1&FR5)	scratch on cut egg white (FR2)	used tool surface aft. 3min. eggshell cutting (FR3)
Grater height	1	anodized aluminum	30	0.2	dent first	1	0.095 [mm³/s]	visible	no wear
	2		20	0.2	dent first	2	0.089	visible	no wear
	3		20	0.2	bulge first	3	0.074	not visible	no wear
	4		5	0.2	bulge first	4	0.010	not visible	no wear
	5		5	0.1	bulge first	5	0.025	slightly visible	no wear
	6	steel*1	20	0.2	bulge first	6	0.069	not visible	local wear on bulges
	7	brass*2	20	0.2	dent first	7	0.008	visible	hard wear on bulges
	8	brass	-20	0.2	dent only	8	0	slightly visible	-
	9	brass	0	-	flat	9	0	not visible	-
Diamond height	10	nickel plated brass	30	-	-	10	0.15	visible	-
	11		20	-	-	11	0.13	visible	-
	12		15	-	-	12	0.10	visible	-
	13		10	-	used in surg.	13	0.087	not visible	particles dropped
	14		7	-	-	14	0.072	not visible	particles dropped

Anodized aluminum tool has no wear.

Surfaces of cut workpieces



Calculated gap heights ($H_{min}R$) by Herrebrugh equations.

tool R=1mm, 200,000rpm E=100GPa, v=0.2 down force w=0.2N, L=2mm	$H_{min}R^*$ at line contact	$H_{min}R^*$ at point contact
(a) eggshell water E=10GPa v=0.2	0.065µm	0.019µm
(b) egg white water E=20kPa v=0.4	10µm	5.0µm
(c) egg white air E=20kPa v=0.4	0.93µm	0.37µm

*1 $H_{min} = 3.10 U^{0.6} W^{-0.2}$ $U = \eta u / (ER)$, $W = w / (ERL)$
*2 $H_{min} = 2.80 U^{0.65} W^{-0.21}$ $U = \eta u / (ER)$, $W = w / (ER^2)$
 $\eta = 0.001 Pa \cdot sec$ (water), $0.000018 Pa \cdot sec$ (air)

Young's ratio : 300kPa of a live pig blood vessel, 20 to 50kPa of boiled egg white.

Simulated pressures on grater tools

